



SELECTED OPPORTUNITIES IN HAEMATOLOGY

Free α -hemoglobin as biomarker for diagnosing and stratifying severity of β -thalassemia and other hemoglobin-related disorders (BIO08398)

FREE α -HEMOGLOBIN AS BIOMARKER FOR DIAGNOSING AND STRATIFYING SEVERITY OF β -THALASSEMIA AND OTHER HEMOGLOBIN-RELATED DISORDERS (BIO08398)

Stage:
Pre-Analytic
Validation

Product factsheet

▶ Biomarker:

- ◆ Free α -hemoglobin (α -Hb)

▶ Technology:

- ◆ ELISA, HTRF

▶ Sample:

- ◆ Blood

▶ Information:

- ◆ Diagnostic
- ◆ Patient Stratification
- ◆ Prognosis

▶ Scientific and Clinical Rationale:

- ◆ The normal development of red blood cells requires a coordinated synthesis of the hemoglobin (Hb) subunits, the α - and β -globins in the case of adult hemoglobin (Hb A). Unlike the β -hemoglobin chains (β -Hb) which are soluble and form homologous tetramers, the free α -hemoglobin chains (α -Hb) are highly instable, and when in excess, form precipitates and act as active oxidants causing apoptosis and inefficient erythropoiesis.
- ◆ The molecular chaperone of α -Hb, the « Alpha-Hemoglobin Stabilizing Protein » (AHSP) specifically binds to α -Hb to form a stable soluble heterodimer but not to the β -Hb or to tetrameric Hb A.

▶ POC:

- ◆ Ctrl healthy donor (n=50); vs β -thalassemic (n=58) or α -thalassemic (n=24) patients
- ◆ Free α -Hb levels are higher in β -thalassemic patients (n=20) compared to α -thalassemic patients (n=6) or healthy subjects without Hb anomalies (n=28).
- ◆ Additional prospective cohort ALPHAPOOL: β -thalassemic (n=58) or α -thalassemic patients (n=24) or healthy subjects (=50)

FREE α -HEMOGLOBIN AS BIOMARKER FOR DIAGNOSING AND STRATIFYING SEVERITY OF β -THALASSEMIA AND OTHER HEMOGLOBIN-RELATED DISORDERS (BIO08398)

Product factsheet

Stage:
Pre-Analytic
Validation

Clinical State and Market Opportunity

▶ Clinical State:

- ◆ Epidemiology:
 - ◆ Affects a prevalent population of around 440000 patients worldwide¹
- ◆ Treatment:
 - ◆ Blood transfusion, Iron Chelation therapy, hydroxycarbamide...

▶ Opportunity:

- ◆ Variety of other potential indications:
 - ◆ Tracking extranumerous α genes: various duplicated or triplicated genotypes
 - ◆ Suspicion of α^0 -thal trait in association with β -thalassaemia trait
 - ◆ Normocytic cobalamin or folate deficiency (masking thalassaemia trait)
 - ◆ Monitoring therapies (inducers of the synthesis of γ -globin, Hydroxycarbamide, Sotatercept...)
 - ◆ Diagnostic orientation needed as an emergency in genetic counselling
 - ◆ Identification of the risk of major β -thalassaemia disease during pregnancy and the diagnosis of difficult symptomatic γ -thalassaemias of newborn

Unique Selling Points

▶ Priority or Patent:

- ◆ EP09 305 352.8 on 2009/04/24
- ◆ PCT/EP2010/055479 on 2010/04/23
- ◆ Granted in US, EP(FR, ES, DE, GB, IT)

▶ Scientific Publication(s):

- ◆ Am J Hematol., 2017 Oct, *Vasseur C. et al.*, doi: 10.1002/ajh.24835
- ◆ BJH, 2017 Apr, *Vasseur C. et al.*, doi: 10.1111/bjh.14800
- ◆ Am J Hematol., 2011 Feb, *Vasseur C. et al.*, doi: 10.1002/ajh.21918

▶ Development opportunities

- ◆ Ongoing research open for partnering


¹ Worldwide: 2015 prevalence of thalassemias; Lancet; V.388; No.10053; 10/8/16; DOI:10.1016/S0140-6736(16)31678-6

FREE α -HEMOGLOBIN AS BIOMARKER FOR DIAGNOSING AND STRATIFYING SEVERITY OF β -THALASSEMIA AND OTHER HEMOGLOBIN-RELATED DISORDERS (BIO08398)

Proof of concept

Biomarker detection: Use of GST-tagged AHSP for specific measurement of free α -Hb

- After expression and purification, 400 μ g of GST-AHSP was bound to the microspin column containing 50 μ L of Glutathione Sepharose 4B and was then incubated with hemolysates of patients to capture the free α -Hb. After a 30 minute incubation at 4°C under gentle agitation, the microspin column was washed five times with PBS (150 mM NaCl, 10 mM Na₂HPO₄, pH 7.4) and the bound proteins were eluted by 200 μ L glutathione buffer (10 mM reduced glutathione in 50 mM Tris-HCl buffer at pH 8.0). The α -Hb containing in the elution fraction was quantified by spectrophotometry at 414 nm ($\epsilon = 125 \text{ mM}^{-1}\cdot\text{cm}^{-1}$) with a HP 8453 spectrophotometer. The total quantity of subunits of Hb in 1 mL of hemolysate was also determined at 414 nm after a 400 fold dilution. The fraction of free α -subunits is thus simply the ratio of absorption of eluted α -subunits from the column on the absorption of total subunits of Hb and is reported in ng/mg equivalent to ppm.




FREE A-HEMOGLOBIN AS BIOMARKER FOR DIAGNOSING AND STRATIFYING SEVERITY OF β -THALASSEMIA AND OTHER HEMOGLOBIN-RELATED DISORDERS (BIO08398)

Proof of concept

► **Biomarker Discovery:** α -Hb levels are higher in β -thalassemic patients compared to α -thalassemic patients or healthy subjects without Hb anomalies.

- ◆ Using the described method, levels of α -Hb relative to total Hb A input are measured
- ◆ β -thalassemic patients (n=20) compared to α -thalassemic patients (n=6) or healthy subjects without Hb anomalies (n=28)





FREE A-HEMOGLOBIN AS BIOMARKER FOR DIAGNOSING AND STRATIFYING SEVERITY OF B-THALASSEMIA AND OTHER HEMOGLOBIN-RELATED DISORDERS (BIO08398)

Proof of concept

► **Pre-Analytic Validation: A prospective study in α -thalassemic and β -thalassemic patients.**

- ◆ Comparison of α -Hb pool between controls (n = 38) and patients with β -thal (n = 57) and α -thal (n = 23).
- ◆ (A) The α -Hb dosing applied to β - and α -thal patients compared with controls without Hb abnormality.
- ◆ (B) Enlargement of the area between 10 and 150 ppm.
- ◆ (C) Sensitivity, specificity and accuracy of the test for predicting individual carrying β -thal comparing of using direct sequencing as gold standard. α -Hb pool is expressed in ppm, equivalent to ng α -Hb/mg of total Hb subunits per ml of haemolysate.



(C)

α -Hb pool	Direct sequencing of genes for β -thal		Total
	Positive	Negative	
ppm > 120	49	0	49
ppm < 120	8	61	69
Total	57	61	118
Sensitivity	86% (49/57)		
Specificity	100% (61/61)		
Accuracy	93% (110/118)		